

Product description

for AffiSep® Glycoprotein standards*

Product description:

Glycoprotein standards consist of qualified glycoproteins to be dissolved in adsorption buffers for the monitoring of lectin column performance and activity. The standards should be used before the application of any samples to the column and after 5 to 10 runs or after a long period of storage in order to determine any loss of activity.

Cat. No.	GP-Standard	Amount
071010	Asialofetuin	5 mg
071011		10 mg
071020	Fetuin	5 mg
071021		10 mg
071040	Glucose Oxidase	5 mg
071041		10 mg
071050	Porcine Thyroglobulin	5 mg
071051		10 mg
071060	Ovalbumin	5 mg
071061		10 mg
071070	Transferrin	5 mg
071071		10 mg
071080	α1-Glycoprotein	5 mg
071081		10 mg

Storage:

- Standards should be stored at 4 °C.

* for research, laboratory and in vitro use only; not for drug or diagnostic use, food or food additives, household or other uses.

Please contact for orders and technical support:

GALAB Technologies GmbH
Tel: +49 / (0)4152 / 889 400
Fax: +49 / (0)4152 / 889 401

Date: 25.08.2009

info@galab.de / www.galab.de